

BONE MARROW BIOPSY & ASPIRATION



Your partner in Biopsy and Special Needles

MED-B

BONE MARROW BIOPSY AND ASPIRATION SYSTEM



- Ergonomic handle
- Luer Lock connection for syringe
- Blunt probe with cm markings to verify specimen length within the needle cannula and to facilitate specimen expulsion
- Available in explant/transplant version
- 5 year shelf life

GAUGE		LENGTH (MM)	
8		08 100	08 150
11		11 100	11 150
13	13 070*	13 100*	13 150*

Other lengths and gauge sizes available as special orders; minimum quantities may apply.
*13 ga "P" Versions do not include an extraction cannula.

MED-L

BONE MARROW BIOPSY AND ASPIRATION SYSTEM
WITH SPECIAL EXTRACTION CANNULA



- Ergonomic handle
- Separate Extraction Cannula to harvest the specimen with minimal crush artifact
- Luer Lock connection for syringe
- Blunt probe with cm markings to facilitate specimen expulsion and to verify specimen within the needle cannula
- 5 year shelf life

GAUGE		LENGTH (MM)	
8		08 100	08 150
11		11 100	11 150
13	13 070*	13 100*	13 150*

Other lengths and gauge sizes available as special orders; minimum quantities may apply.
*13 ga "P" Versions do not include an extraction cannula.

MED-J
BONE MARROW BIOPSY AND ASPIRATION SYSTEM

- Luer Lock connection for syringe
- Blunt probe with cm markings to verify specimen length within the needle cannula and to facilitate specimen expulsion
- 5 year shelf life



GAUGE	LENGTH (MM)	
8	08 100	
11	11 100	11 150
13	13 050	13 090

Other lengths and gauge sizes available as special orders; minimum quantities may apply.

MED-I
BONE MARROW ASPIRATION SYSTEM

- Ergonomic handle for a comfortable procedure and a firm grip
- Triple sharpened to allow an atraumatic, fast and easy penetration into the bone marrow cavity
- Adjustable depth control
- Luer Lock connection for syringe
- Available in explant/transplant version
- 5 year shelf life



GAUGE	LENGTH (MM)					
14	14 030	14 040	14 050	14 070	14 090	14 110
15	15 030	15 040	15 050	15 070	15 090	15 110
16	16 030	16 040	16 050	16 070	16 090	16 110
18	18 030	18 040	18 050	18 070	18 090	18 110

Other lengths and gauge sizes available as special orders; minimum quantities may apply.

ORDERING INFORMATION

MED-B

MED-B DIAMOND TIP

X	B	GAUGE	LENGTH	0	0

Single box: 10 PCS.

MED-L

MED-L “S” VERSION DIAMOND TIP

X	L	GAUGE	LENGTH	0	S

MED-L “P” VERSION FISH MOUTH TIP

X	L	GAUGE	LENGTH	0	P

Single box: 10 PCS.

MED-J

MED-J

X	J	GAUGE	LENGTH	0	0

Single box: 10 PCS.

MED-I

MED-I ASPIRATION TIP

X	I	GAUGE	LENGTH	0	1

Single box: 10 PCS.



MEDAX S.R.L. Unipersonale

Via Sandro Pertini, 4 • 41039 • San Possidonio (MO) • Italy

Company direct No. : +39 0535 1812757 • customercare@medax.it • www.medax.it • www.medax-usa.com

Vat N. /Fiscal Code N. Iscriz. Reg. Impr. MN 02669860369 • N. REA: MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates.
All products undergo intensive clinical testing and are fully EC and FDA approved.