

PROSTATE BIOPSY



Your partner in Biopsy and Special Needles

MEDONE

AUTOMATIC DISPOSABLE BIOPSY INSTRUMENT



- One hand usage
- 22mm penetration depth, 19mm sample notch
- Dual firing triggers in proximal and distal positions
- Echogenic tip for ultrasound guided procedures
- Centimeter markings

GAUGE	COLOUR	LENGTH (MM)		
18	Pink ●	18 200	18 250	18 300

Ordering information:

MEDONE

Single box: 10 PCS.

M

O

GAUGE

LENGTH

0

0

MEDCUT

AUTOMATIC DISPOSABLE BIOPSY INSTRUMENT

- 22mm penetration depth, 19mm sample notch
- Firing trigger in proximal position
- Extra sharp tip for safe insertion and penetration
- Echogenic tip for ultrasound guided procedures
- Centimeter markings



GAUGE	COLOUR	LENGTH (MM)		
18	● Pink	18 200	18 250	18 300

Ordering information:

MEDCUT

Single box: 10 PCS.

M

C

GAUGE

LENGTH

0

0

CHIBA

FINE NEEDLE ASPIRATION

Fine needle for cytological aspiration biopsy or local anesthesia administration

- Transparent Luer Lock hub
- Sharpened cannula for safe and painless insertion
- Internal and external echogenic treatment for correct positioning under ultrasound guidance
- Centimeter markings

GAUGE	COLOUR	LENGTH (MM)			
22	● Black	22 180	22 200	22 250	22 300



Ordering information:

CHIBA NEEDLE WITH ECHOGENIC TREATMENT



Single box: 20 PCS.

MED-GUIDE 1

DISPOSABLE ENDOCAVITY BIOPSY NEEDLE GUIDE



Disposable endocavity biopsy needle guide used in ultrasound guided procedures such as transrectal prostate biopsies, transvaginal biopsies and transvaginal oocyte recovery. Compatible with: GE Healthcare E8C, E8C-RS, E8CS, E72, E721, IC5-9-D, IC5-9H, MTZ, MZ, ECM, E610P and Samsung VR5-9 transducers

- Precise needle placement
- Smooth profile ensures patient comfort
- Requires no disinfecting
- Compatible with needles 16G or thinner

Ordering information:

MED-GUIDE 1



Single box: 24 PCS.

UNIVERSAL

BIOPSY NEEDLE USABLE WITH MANY
COMMONLY USED REUSABLE HANDLE BIOPSY SYSTEMS

Examples are: Pro-Mag™ 2.2 Handle
Biopty® Handle
ABS Handle
MBD-23 Handle

- Echogenic treatment for ultrasound guided procedures
- Centimeter markings



UNIVERSAL DUO

BIOPSY NEEDLE USABLE WITH PRO-MAG™ ULTRA 2.5 &
COOK® Bx™ REUSABLE HANDLE BIOPSY SYSTEMS


- Echogenic treatment for ultrasound guided procedures
- Centimeter markings



UNIVERSAL PLUS

BIOPSY NEEDLE USABLE WITH MAGNUM®
REUSABLE HANDLE BIOPSY SYSTEM

- Echogenic treatment for ultrasound guided procedures
- Centimeter markings



UNIVERSAL / UNIVERSAL DUO / UNIVERSAL PLUS

GAUGE	COLOR CODE	LENGTH (MM)		
18	● Pink	18 200	18 250	18 300

Ordering information:

UNIVERSAL

U	N	GAUGE	LENGTH			0	0
---	---	-------	--------	--	--	---	---

UNIVERSAL DUO

U	D	GAUGE	LENGTH			0	P
---	---	-------	--------	--	--	---	---

UNIVERSAL PLUS

U	P	GAUGE	LENGTH			0	0
---	---	-------	--------	--	--	---	---

Single Box: 10 PCS

**MEDAX S.R.L. Unipersonale**

Via Sandro Pertini, 4 • 41039 • San Possidonio (MO) • Italy

Company direct No. : +39 0535 1812757 • customercare@medax.it • www.medax.it • www.medax-usa.com

Vat N. /Fiscal Code N. Iscriz. Reg. Impr. MN 02669860369 • N. REA: MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates.
All products undergo intensive clinical testing and are fully EC and FDA approved.